

Product name:	FR3 Battery Charger	
Identification code(s)	861394	
Total weight (in Kg)	0.44 kg -0.52 kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	22100 Bothell Everett Highway
	Zip code:	Bothell, WA 98021
	Country:	United States






Recycle Info	Items:	Location
Special attention 	Remove rubber feet (x4) on the back of charger.	Figure 1
	Remove screws on the back of charger to open.	Figure 2
	Contains circuit board.	Figure 3
 Fluids / Gases	Items:	Location
	N/A	
Batteries  To be Removed	Type:	Location
	N/A	
Hazardous  To be Removed	Substances:	Location
	N/A	


FIGURE 1



Rubber Feet (4X)



FIGURE 2



SCREWS



FIGURE 3

SCREWS

PCA

