

The Achieva 1.5T – Diamond Select helps you to do more in every way

Faster exams, a wide application range and enhanced workflow are just some of the results. Exceptional quality in challenging applications, across all patients.

A smart investment by any measure. The Achieva 1.5T – Diamond Select provides high patient friendliness, functionality and versatility. This high field MR scanner is fully customizable to fit your working style and situation.

Key advantages

- Drive clinical performance with excellent magnet homogeneity, Achieva gradients and Examcards
- Enhanced patient experience with Spacious design, SofTone and (optional) ComforTone
- Ensured economic value wtih scalable platform for now and in the future



Drive clinical performance

- Excellent magnet homogeneity and the Achieva 1.5T gradient systems enable the performance for today's demanding and advanced clinical applications
- ExamCards make it easy to download and use protocols and therefore shorten exam time and reduce training requirements

Enhance patient experience

- Appealing spacious design that provides high patient access
- SofTone technology and ComforTone (optional) reduce patient-perceived acquistic noise

Ensure economic value for now and in the future

- Scalable platform from which you choose your gradients and channel amount for optimal performance
- The advanced architecture allows you to add clinical packages as your needs change and to expand your application mix and referral base when required.

The Achieva 1.5T – Diamond Select provides an enhanced growth path by clinical strength, patient comfort and economical value for now and the future.

The Achieva 1.5T is available with Pulsar gradient (amplitude 33 mT/m, slew rate 80 mT/m/ms) Nova gradient (33/160) and Nova Dual gradient (33/160 or 66/80). The Nova Dual gradient system delivers the advanced performance levels required for demanding clinical applications.

First-rate care

Philips partners with you to help you get the most out of your investment. Our global network of skilled personnel and worldwide spare parts availability helps maintain high uptime that supports you in delivering reliable imaging services to your patients. Philips has been recognized by Best in KLAS' for the last several years as a performance leader.

2013 and 2014 Best in KLAS Award: Medical Imaging & Infrastructure report can be found at: www.klasresearch.com

