



PHILIPS
RESPIRONICS

DreamStation Go

Portable PAP
therapy system

Give your patients convenience **without compromise**

Powerful, portable, performance

At half the size of our System One and DreamStation platforms, DreamStation Go offers the same therapy and comfort features used by over 5 million CPAP users.

Patient preferred tubing¹

Our 12mm micro-flexible tubing is designed to provide your patients with added freedom of movement while they sleep. Engineered to minimize bulk and improve mobility, our 12mm tubing is 58% smaller than our 22mm performance tubing. It offers enhanced flexibility, portability, and connectivity. Click-connections on each end of the micro-flexible tubing makes connecting and disconnecting the device and mask simple.

Innovation for patients on the go

DreamStation Go PAP provides the portability and packability travelers want and the effective therapy they need.

Proven Flex technology

Our clinically-proven Flex technologies enhance traditional CPAP therapy and are an excellent way to improve patient comfort. DreamStation Go utilizes the same Flex technology as our DreamStation and SytemOne PAP devices.

DreamMapper cross-device connectivity

Users now have the ability to track and monitor their combined progress across both home and portable PAP devices in one convenient dashboard on their smartphone, tablet or computer.

Powerfully convenient features

- Intuitive color touchscreen offers simple setup and easy navigation
- Built-in power supply reduces parts managed by patients
- USB port for easy charging of mobile devices
- Bluetooth is standard with every device
- Complies with FAA requirements for in-flight use²

33% fewer components³

The DreamStation Go requires 33% fewer components to pack than its closest competitors.

¹Head-to-Head trial of DreamStation 15mm tubing vs. DreamStation 12mm micro-flexible tubing, 2016.

²RTCA/DO-160G section 21, category M.

³When compared to HDM Z1 and Somnetics Transcend devices, (2016).

DreamStation Go Auto specifications

Pressure range	4 to 20 cm H ₂ O
Pressure mode(s)	CPAP, APAP
Flex pressure relief	1 to 3
Flex modes	C-Flex, C-Flex+, A-Flex
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable, Smart Ramp
Dimensions	15 cm L x 15 cm W x 5.8 cm H
Weight	844 g (1.86 lb.)
Filters	Reusable pollen, combination fine
Device set-up	Touchscreen color display
Data storage capacity	On-board: > 1 year MicroSD card (optional): lifetime
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100-240 VAC, 50/60 Hz
Warranty	Two years

Ordering information

Flow Generators	Part number
DreamStation Go Auto CPAP	GBG502S15
Power	Part number
DreamStation Go overnight battery	1133281
Travel kits	Part number
DreamStation Go small travel kit ¹	1133275
DreamStation Go medium travel kit ²	1133276
Consumables/Replacement	Part number
DreamStation micro-flexible 12mm tubing	PR12
Power cord 6ft white	1132551
Power cord 10ft white	1132552
Beauty panel - front	1133741
Beauty panel - rear	1133742
Reusable pollen filter 2-pk	1133743
Disposable 30-day fine filter 6-pk	1133744

Accessories create opportunities for you



Small travel kit

Kit includes an easy-to-pack protective case for the DreamStation Go and a separate, washable bag for the mask and tube.



Medium travel kit

Kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes separate, washable bag for the mask and tube.



Air filters

Your choice of a reusable pollen or disposable 30-day fine filter.



Overnight battery

Compact and lightweight with a recharge time of just 4-5 hours.

¹Small travel kit includes an easy-to-pack protective case for the DreamStation Go and a separate washable mask and tube bag.

²Medium travel kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes a separate washable mask and tube bag.

