

Mask and accessories materials list

Includes all lab versions except where noted.

No parts contain natural rubber latex.

Product	Material
Amara	
Frame/hub	Polycarbonate
Cushion	Silicone
Forehead pad	Silicone
Exhalation elbow	Polycarbonate
Exhalation chip	Acetal copolymer
Elbow swivel	Polycarbonate
Split washer	Polycarbonate
Forehead adjuster/bracket	Polycarbonate
Headgear clip	Thermoplastic polyester elastomer
Headgear	UBL/Urethane foam/nylon and Lycra

TrueBlue	
Faceplate	Polycarbonate
Gel cushion	Polyurethane gel and polyester film
Gel forehead pad	Polyurethane gel and polyester film
Silicone spring ring	Polycarbonate
Silicone spring	Silicone
Silicone spring male ring	Polyether Imide
Sealing flap	Silicone
Sealing flap ring	Polycarbonate
Exhalation elbow	Polycarbonate
Exhalation chip	Acetal copolymer
Exhalation swivel	Polycarbonate
Split washer	Polycarbonate
Forehead bracket	Thermoplastic polyester elastomer
Talon clip	Thermoplastic polyester elastomer
Headgear	UBL/Urethane foam/nylon and Lycra

GoLife	
Frame	Thermoplastic polyester elastomer
Cushion	Silicone
Tubing	Ethylene octane copolymer
Exhalation elbow and swivel	Polycarbonate
Headgear and frame padding	UBL/Urethane foam/Lycra

Product	Material
ComfortGel Blue	
Faceplate	Polycarbonate
Gel cushion	Urethane gel/EVA
Gel cushion outer shell	Urethane film
Outer flap	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Forehead support arm and bracket	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear (premium)	UBL/Urethane foam/Lycra

ComfortGel	
Faceplate	Polycarbonate
Gel cushion	Urethane gel
Gel cushion outer shell	Urethane film
Outer flap	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Forehead support arm and bracket	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear (premium)	UBL/Urethane foam/nylon and Lycra

EasyLife	
Faceplate	Polycarbonate
Inner seal cushion	Silicone
Outer support cushion	Silicone
Diaphragm cushion	Silicone
Exhalation elbow	Polycarbonate
Standard elbow (SE)	Polycarbonate
Elbow split washer	Polycarbonate
Headgear	UBL/Urethane foam/Lycra

Product	Material
ComfortGel Full	
Faceplate	Polycarbonate
Gel cushion	Urethane gel
Gel cushion outer shell	Urethane
Outer flap (SST)	Silicone
Cushion flap ring	Polycarbonate
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Silicone
Entrainment valve flapper	Silicone
Forehead support arm and bracket	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear (premium)	UBL/Urethane foam/Lycra

ComfortGel Blue Full	
Faceplate	Polycarbonate
Gel cushion	Urethane gel
Gel cushion outer shell	Urethane
Outer flap (SST)	Silicone
Cushion flap ring	Polycarbonate
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Exhalation chip	Acetal copolymer
Elbow split washer	Polycarbonate
Entrainment valve flapper	Silicone
Forehead support arm and bracket	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear (premium)	UBL/Urethane foam/Lycra

FullLife	
Frame	Thermoplastic polyester elastomer
Outer flap (SST)	Silicone
Entrainment elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Entrainment valve flapper	Silicone
Headgear	UBL/Urethane foam/Lycra and nylon

FitLife	
Faceplate	Polycarbonate
Cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Faceplate valve hub	Polycarbonate
Accessory swivel	Polycarbonate
Elbow split washer	Polycarbonate
Entrainment valve flapper	Silicone
Port cap	Polycarbonate
Headgear clips	Acetal
Headgear	UBL/Urethane foam/Lycra
Headgear (lab version)	Laminate composition of nylon and polyester polyurethane foam

Product	Material
OptiLife	
Frame	Thermoplastic polyester elastomer
Cushion (pillows and CradleCushion)	Silicone
Tubing	Ethylene vinyl acetate (EVA)
Swivel	Polycarbonate
Headgear	UBL/Urethane foam/Lycra

ComfortFusion	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate and/or PETG
Forehead support arm and bracket	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear	UBL/Urethane foam/Lycra

ComfortLite 2	
Cushion (pillows, Simple Cushion, Direct Seal)	Silicone
Simple Cushion elbow	Polycarbonate
Direct Seal cushion cradle	Polycarbonate
Direct Seal cushion support band	Aluminum
Crown swivel	Polycarbonate
Flex tube	EVA copolymer
Pivot tube insert	Polycarbonate
Vertical adjustment nut	Polypropylene
Horizontal adjustment tube	Polypropylene
Pivot mount strap	Polycarbonate
Pivot tube	Polypropylene
Headgear halo	Nylon/polyester
Headgear padding	Nylon/Lycra laminated foam
Headgear strap	Nylon/Lycra laminated foam

ComfortFull 2	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Entrainment valve flapper	Silicone
Port cap	PVC
Headgear (premium)	UBL/Urethane foam/Lycra
Headgear clips	Thermoplastic polyester elastomer

Product	Material
ComfortSelect	
Faceplate	Polycarbonate
Cushion	Silicone
Forehead cushion	Silicone
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate
Headgear clips	Thermoplastic polyester elastomer
Headgear (deluxe)	Velstretch/Lycra laminate

ComfortClassic	
Faceplate	Polycarbonate
Cushion	Silicone
Gel forehead pad	Urethane gel encased in a polypropylene sheet
Exhalation elbow and swivel	Polycarbonate
Elbow split washer	Polycarbonate and/or PETG
Port cap	Flexible vinyl (PVC)
Headgear (deluxe)	Velstretch/Lycra laminate

Respironics Profile Lite, Profile Lite Youth, Small Child Profile Lite	
Faceplate	Polycarbonate
Cushion outer shell	Urethane
Port cap	Flexible vinyl (PVC)
Exhalation elbow and swivel	Polycarbonate
Spacer	White polyethyl foam
Headgear (deluxe)	Velstretch/Lycra laminate
Child SoftCap headgear	Velstretch/Lycra laminate
Headgear (RS)	UBL/Urethane foam/Lycra

Respironics Simplicity nasal mask	
Faceplate	Polycarbonate
Cushion	Silicone
Tubing	Ethylene vinyl acetate (EVA)
Deluxe softcap	Velstretch/Lycra nylon
Headgear	Neoprene/nylon

Product	Material
Respironics Total face mask	
Faceplate	Polycarbonate
Cushion	Silicone
Seal ring	Thermoplastic polyester elastomer
Valve body	Polycarbonate
Headgear	Neoprene/nylon

Lab masks (LabClassic, LabFull, LabGel, LabSelect)

Materials as listed above for reusable masks, except:

Foam forehead pad	Polyurethane foam
Disposable headgear	Polyurethane foam lined with nylon cloth (40 denier tricot)

Headgear

Headgear (premium)	UBL/Urethane foam/Lycra
Headgear (RS)	UBL/Urethane foam/Lycra
Headgear (deluxe)	Velstretch/Lycra laminate
SimpleStrap:	
Mesh	Polyester
Strap loop/hook	Velstretch/nylon hook
Softcap:	
Cap	Lycra
Strap loop/hook	Velstretch/nylon hook

Tubing

Standard gray tubing:	
Cuffs	Silicone rubber or equivalent thermoplastic elastomer
Wall	Polyester elastomer or ethylene octene copolymer
Helix	Polyester elastomer or LDPE or HDPE or equivalent
Performance tubing:	
Cuffs	Thermoplastic elastomer
Wall	Ethylene octene copolymer
Helix	Polypropylene

Respironics and the names of Respironics products included here are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. Names of non-Respironics products are trademarks of their respective manufacturers.

Please visit www.philips.com/respironics



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

edoc JJ 06/18/12 MCI 4104975

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics