

A clear advantage in ERCP procedures

Philips MultiDiagnost Eleva system

For ERCP (Endoscopic Retrograde Cholangio Pancreatography) procedures it is essential to obtain a detailed visualization of complex vasculature from multiple viewing angles. The advanced high resolution from the MultiDiagnost Eleva allow the visualization of small anatomical details, while 3D-RX imaging provides better insight in the complex anatomy of the biliary system. The MultiDiagnost Eleva system with the ERCP package provides all of this and more.

PHILIPS
sense **and** simplicity

ERCP package

The special package for ERCP examinations is designed to perform the most complete ERCP procedures and includes 3D-RX, which provides excellent image quality attained with 180° rotation. This package can also be used for GI (Gastrointestinal) diagnostic and interventional examinations.

Provides high resolution images for better results

- High resolution 2K images visualize small anatomical details and support ERCP interventions making it easy to get the desired result.
- The projection flexibility offered by the C-arc of the MultiDiagnost Eleva supports attaining excellent anatomical views from various angles.

Enhances patient and staff comfort during ERCP procedures

- The C-arc moves around the table, so that the patient can remain stationary during lengthy interventional procedures, increasing comfort for the patient and gastroenterologist.
- For optimal access to the patient during an ERCP procedure, a left suspended table is available. This allows easy access for anesthesia needs and interventional treatment.

Enhances insight during ERCP procedures

- 3D-RX reconstructions provide enhanced insight in the complex anatomy of the biliary system and enable navigation.
- 3D-RX images are available within seconds during treatment.

Features at a glance


- Superb high resolution images for enhanced results
- C-arm movement for patient comfort
- Projection flexibility for optimal treatment view
- Left suspension of table for optimal patient access
- 3D images for confidence in navigation
- DoseWise for dose management



ERCP: 3D Biliary tree



ERCP: 3D Pancreas



ERCP: Pancreas cyst

Please visit www.philips.com/multidiagnosteleva



© 2012 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 87921 * SEP 2012