


Unique as every face

Respironics Profile Lite adapts to facial features for a custom fit every time.

This fully-customizable nasal gel mask fits comfortably right off the shelf, yet is designed to be easily contoured to the unique shape of any wearer's face in just minutes. Thanks to a soft, light cushion and Respironics' proven gel technology, the Profile Lite nasal gel mask offers an unmatched, personalized seal. A simple customization process allows you to quickly mold the mask to a wearer's facial contours. Profile Lite provides the simplest, most sensible way to achieve the highest level of comfort and performance for all wearers—men, women, and even your pediatric patients. It's just one more way we're helping you fit 100% of your patients, 100% of the time.





Swivel allows tubing to move freely.

Built-in silent exhalation port disperses the air for quiet operation.

Pressure ports allow pressure to be measured at the mask.

Inner layer, when customized, molds to facial contours.



Before molding: standard contour shape



After molding: customized contour shape

Custom fitting Profile Lite is easy. Simply boil, cool, and press the mask into place on your patient's face.

Item				
Profile Lite mask	With headgear	HCPCS code	Mask only	HCPCS code
Petite	1004086	A7034,A7035	1002849	A7034
Small	1004087	A7034,A7035	1002850	A7034
Medium small	1004088	A7034,A7035	1002370	A7034
Medium	1004089	A7034,A7035	1002371	A7034
Medium wide	1004110	A7034,A7035	1002372	A7034
Large	1004111	A7034,A7035	1002373	A7034
Large narrow	1004112	A7034,A7035	1002374	A7034
Profile Lite Youth Mask with standard deluxe headgear	1047339	A7034,A7035	N/A	N/A
Parts and accessories			Part number	HCPCS code
Deluxe headgear			1002800	A7035
Deluxe headgear, small			1018755	A7035
Exhalation elbow kit			1004395	A7045
Gel retrofit forehead spacer, small			1004292	N/A
Gel retrofit forehead spacer, large			1004291	N/A
Gel retrofit forehead spacer, extra large			1008869	N/A
Gel retrofit forehead spacer, assortment pack (3 small, 1 large, 3 extra large)			1009857	N/A

Respironics and Respironics Profile Lite Nasal Gel Mask are trademarks of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Philips Healthcare is part of Royal Philips Electronics

Europe, Middle East, Africa: +33 1 47 52 30 00

Asia Pacific: +852 3194 2280 Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 +1 724 387 4000

+1 800 345 6443

Hoech JJ 2/17/09 MCI 4102139 PN 1059961