



Your oral/enteral solution for the NICU

Oral/enteral syringe, sterile

In response to the Joint Commission's Sentinel Event Alert, Issue 36, April 2006 – Tubing misconnections, a persistent and potentially deadly occurrence – Philips Children's Medical Ventures has collaborated with healthcare professionals to develop a line of oral/enteral syringes for use with neonatal, infant and pediatric patients throughout the hospital.

Key advantages

- Sterile and individually packaged
- Bright orange coloring to alert caregivers to enteral only use
- Can help avoid dangerous feeding and medication errors that can pose a significant risk to the infant and the hospital

PHILIPS

Children's  Medical Ventures
sense and simplicity

Oral/enteral syringes

The oral/enteral solution for neonatal, infant and pediatric patients

Philips Children's Medical Ventures' complete line of oral/enteral syringes has been developed to help reduce feeding errors and to avoid tubing misconnections that can pose a significant risk to the patient and the hospital.

Oral/enteral syringes come in seven different sizes (1 ml, 3 ml, 5 ml, 10 ml, 20 ml, 30 ml and 60 ml) with additional features that include highly visible labeling that reads "ORAL/ENTERAL ONLY," a bright orange plunger that immediately distinguishes it from medication syringes, and a syringe tip that will not securely attach to a standard luer lock connector. In addition, the 1 ml syringe – the size commonly used to administer medication in the hospital setting – is a highly noticeable orange color, with labeling that reads "ORAL/ENTERAL ONLY."

Features

- Sterile, individual packaging and barcode labeling
- Bright orange plunger alerts clinician that syringe is for oral use only
- Large, highly visible labeling on side of barrel to easily distinguish enteral feeding, even during gavage feeding when plunger is removed
- Pointer plunger tip on 1 ml, 3 ml, and 5 ml sizes helps minimize medication or expressed breast milk waste
- One-piece barrel construction assures syringe is more secure than adhesive-secured tips and eliminates separation during filling
- Will not allow connection to a luer lock connector
- Highly visible color of 1 ml oral syringe distinguishes it as an orally administered medication syringe and is ideal for administering Sweet-Ease sucrose solution
- Compatible with the most commonly used hospital syringe pumps



Highly visible labeling on side of barrel marked ORAL/ENTERAL ONLY and bright orange plunger easily distinguishes syringe for feeding.



Syringe caps

- Airtight cap will not allow water in, or breast milk out, during prep, storage and warming
- Cap stands on end for stability when prepping milk on a counter or other surface
- Bright orange color prevents caps from getting lost in bedding



One-piece design eliminates separation during filling and pointed plunger tip on 1 ml, 3 ml, and 5 ml sizes allows for all contents to be delivered.



Individual packaging protects syringe during shipping and storage and is labeled with a barcode for easy inventory management.

Ordering information

Item no.	Description	Quantity
1062677	1 ml oral/enteral syringe, sterile	100/box
1062678	3 ml oral/enteral syringe, sterile	100/box
1062679	5 ml oral/enteral syringe, sterile	100/box
1062680	10 ml oral/enteral syringe, sterile	100/box
1062681	20 ml oral/enteral syringe, sterile	100/box
1062682	30 ml oral/enteral syringe, sterile	100/box
1062683	60 ml oral/enteral syringe, sterile	100/box
1078782	Sample pack, oral/enteral syringe, sterile	1/each

Sweet-Ease is a trademark of Koninklijke Philips Electronics N.V. All rights reserved.

Please visit www.childmed.com



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 71111 * MAR 2011