



# One solution

## Philips ECG cables and leads quick reference guide

One ECG cable solution that is compatible with Philips IntelliVue, CMS and Vital Signs monitors, and HeartStart defibrillator/monitors for hospital-wide standardization.

Philips offers a new cable solution for ECG and respiration monitoring in all patient care areas of your hospital. This allows you to standardize on fewer cables and leads for all your Philips monitors, defibrillators, as well as telemetry, in the ER, OR, ICU, NICU or transport.

Your existing IntelliVue Telemetry 3-, 5- and 6-lead sets are also compatible with the new ECG trunk cables making it simple to convert.

### Features

- Convenient and cost-effective solution
- Available in 3-, 5- and 6-lead including EASI and 10-lead configuration for conventional diagnostic 12-lead ECG
- Lightweight cabling is more comfortable for patient sizes neonates to adults
- New flexible material means fewer tangles for faster application
- Limb lead wires are fully color coded, using AAMI Color Coding to help with rapid identification and placement (according to ANSI/AAMI EC-13 and Code 2 of 60601-2-37)
- Cadmium- and PVC-free for unproblematic waste management

### Patient trunk cable and lead set selection guide: (AAMI Color Coding) 3-lead trunk cable (12-pin connector)

If you have these monitoring devices	Use this trunk cable	Trunk cable length	With one of these lead sets	To replace CMS trunk cables	To replace CMS lead sets
ECG modules: M1001A/B, M1002A/B	M1669A	9 ft. (2.7m)	M1671A (ICU, grabber)	M1500A	M1603A
Multi-Measurement modules: M3000A, M3001A, M3002A X2			M1673A (ICU, snap)	9 ft. (2.7m)	M1605A
Patient monitors: M8102A MP2, M8105A MP5			M1675A (OR, grabber)	(future obsolescence)	M1601A
Vital Signs monitors: 863063 VM4, 863064, 863065 VM6, 863066, 863068 VM8, 863077 VM3			M1622A (mini, grabber)*		M1608A
HeartStart XL defibrillator/monitor: M4735A			M1624A (mini, grabber)*	M1540C	M1609A
HeartStart MRx defibrillator/monitors: M3535A, M3536A				3 ft. (0.9m) (future obsolescence)	
Legacy HP patient monitors: 78352C, 78354C, 78834C					

Note: Third party lead sets are not compatible with the new ECG trunk cables.

\*Miniclips not to be used with defibrillators.

# PHILIPS

**Patient trunk cable and lead set selection guide: (AAMI Color Coding) 5-lead trunk cable (12-pin connector)**

If you have these monitoring devices	Use this trunk cable	Trunk cable length	With one of these lead sets	To replace CMS trunk cable	To replace CMS lead sets
ECG modules: M1001A/B, M1002A/B	M1668A	9 ft. (2.7m)	M1968A (ICU grabber, limb) M1644A (ICU snap, limb) M1973A (OR grabber, limb) M1647A (miniclip, grabber)*	M1520A 9 ft. (2.7m) (future obsolescence)	M1623A M1625A M1621A M1629A
Multi-Measurement modules: M3000A, M3001A, M3002A X2					
Patient monitors: M8102A MP2, M8105A MP5,					
Vital Signs monitors: 863064, 863065 VM6, 863066, 863068 VM8				M1560C 3 ft. (0.9m) (future obsolescence)	
HeartStart XL defibrillator/ monitor: M4735A					
HeartStart MRx defibrillator/ monitors: M3535A, M3536A					
Legacy HP patient monitors: 78352C, 78354C, 78834C					

Note: The M1668A (5-lead trunk cable) can also be used with the new 3-lead sets.

\*Miniclips not to be used with defibrillators.

**Patient trunk cable and lead set selection guide: (AAMI Color Coding) 6-lead trunk cable (12-pin connector)**

If you have these monitoring devices	Use this trunk cable	Trunk cable length	With one of these lead sets	Lead wire length
Patient monitors: M8102A MP2 Opt. C12, M8105A MP5 Opt. C12	M1667A	9 ft. (2.7m)	M1680A (ICU grabber, limb) M1682A (ICU snap, limb) M1684A (OR grabber, limb)	6-lead 5.2 ft. (1.6m) long
IntelliVue Multi-Measurement modules: M3001A Opt. C12 or C18, M3002A X2 Opt. C12 or C18				

Note: The M1667A (6-lead trunk cable) can also be used with the new 3-lead or 5-lead sets.

### Patient trunk cable and lead set selection guide: (AAMI Color Coding) 10-lead trunk (5-lead + 5-lead) 12-pin connector

If you have these monitoring devices	Use this trunk cable	Trunk cable length	With one of these lead sets	To replace IntelliVue trunk cable
IntelliVue Multi-Measurement modules: M3001A Opt. C12 or C18, M3002A X2 Opt. C12 or C18  Patient monitors: M8102A MP2 Opt. C12 M8105A MP5 Opt. C12,  HeartStart MRx defibrillator/monitors: M3535A, M3536A	M1663A*	6.5 ft. (2.0m)	M1968A (ICU grabber, limb) M1976A (ICU grabber, chest) M1644A (ICU snap, limb) M1602A (ICU snap, chest) M1973A (OR grabber, limb) M1979A (OR grabber, chest)	M1949A* 9 ft. (2.7m) (AAMI/IEC)  M1949A trunk cable is not validated for use on HeartStart MRx defibrillator/monitors

Note: For 10-lead applications, choose one 'chest' and one 'limb' lead set.

\*The M1663A and M1949A (10-lead trunk cables) can also be used with the new 3-lead or 5-lead sets.

### Patient trunk cable and lead set selection guide: (AAMI Color Coding) 10-lead trunk (6-lead + 4-lead) 12-pin connector

If you have these monitoring devices	Use this trunk cable	Trunk cable length	With one of these lead sets	Lead wire length
Patient monitors: M8102A MP2 Opt. C12, M8105A MP5 Opt. C12  IntelliVue Multi-Measurement modules: M3001A Opt. C12 or C18, M3002A X2 Opt. C12 or C18	M1665A*	9 ft. (2.7m)	M1680A (ICU grabber, limb) M1532A (ICU grabber, chest) M1682A (ICU snap, limb) M1537A (ICU snap, chest) M1684A (OR grabber, limb) M1557A (OR grabber, chest)	6-lead = limb leads 5.2 ft. (1.6 m)  4-lead = chest leads 3 ft. (1.0m)

Note: For 10-lead applications, choose one 'chest' and one 'limb' lead set.

\*The M1665A (10-lead trunk cable) can also be used with the new 3-lead, 5-lead or 6-lead sets..

### Cable accessories (AAMI Color Coding)

	Model No.	Order No.	With one of these lead sets
<b>Set organizer (shielded)</b>			
3-electrode	M1503A	9898031 03841	M1671A, M1673A, M1675A
4-electrode	M1664A	9898031 50501	M1532A, M1537A, M1557A
5-electrode	M1504A	9898031 03851	M1968A, M1644A, M1973A, M1976A, M1602A, M1979A
6-electrode	M1679A	9898031 50511	M1680A, M1682A, M1684A, M1532A, M1537A, M1557A
<b>Set organizer (unshielded)</b>			
3-electrode	M1636A	9898031 44971	M1622A, M1624A
5-electrode	M1638A	9898031 44981	M1647A

### Other accessories

Description			
Bedsheet clip	M1509A	9898031 03861	
Trunk cable covers (10 per bag)		9898031 48861	For use with M1663A or M1949A 10-lead trunk cables
IntelliVue cable management kit (5 per bag)		9898031 48841	For use to store cables while not connected to patients

**Philips Healthcare is part of  
Royal Philips Electronics**

**Interested?**

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be glad to hear from you.

**On the web**

[www.philips.com/healthcare](http://www.philips.com/healthcare)

**Via email**

[healthcare@philips.com](mailto:healthcare@philips.com)

**By fax**

+31 40 27 64 887

**By mail**

Philips Healthcare  
Global Information Center  
P.O. Box 1286  
5602 BG Eindhoven  
The Netherlands

**By phone**

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +49 7031 463 2254

Latin America

Tel: +55 11 2125 0764

North America

Tel: +1 425 487 7000

800 285 5585 (toll free, US only)



© 2008 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.  
4522 962 35401 \* FEB 2009