

Product name:	Lateral Stand Allura XPER FD10/FD20		
Identification code(s)	4522 901 00140		
Total weight (in Kg)	800 Kg (approx)		
Producer /Manufacturer	Name company:	Philips Medical Systems	
	Address:	Veenpluis 6	
	Zip code:	5684 PC Best	
	Country:	The Netherlands	
	Electronic info:	Lateral Stand_FD10_579.doc	


Recycle info	Substances:	Location
Special attention 	Items:	Location
	Be aware of possible contaminated system parts and materials.	See Figure 1
	For dismantling activities Treatment Facilities must consider the national requirements.	
	For personnel that can come in contact with contaminated material, preventive measurements pursuant to national requirements must taken into account.	
	Removal of units/weights can cause the system to tilt!	
	Removal of units/weights can cause unexpected movements.	
	High voltage parts are marked with	
Hazardous 	Lead (Pb) for soldering	
	Beryllium copper (BeCu)	Some electronics
To be removed		

Note: To facilitate recycling, all plastic parts > 50 grams are marked according to ISO 11469 & ISO 1043


Name: WEEE Lat Stand Allura XPER FD10/FD20		12NC: 4522 901 00140		Version: 00 Doc.Date: 2006-09-01	
Author: W.Bleyen		Pag.: 1 of 3	Group No: 579	Section: -	Status: Approved

Lateral stand

Figure 1




View A



Name: WEEE Lat Stand Allura XPER FD10/FD20		12NC: 4522 901 00140		Version: 00 Doc.Date: 2006-09-01	
Author: W.Bleyen		Pag.: 2 of 3	Group No: 579	Section: -	Status: Approved

View B



Recycling Passport Revision History

Rev	Description	Date
1	First Passport	Sept. 1, 2006

Recycling Passport Approval Signatures

Written By PMS J.Dries

Reviewed By (optional)

Approved By PMS W.Bleyen

Approved By

Name: WEEE Lat Stand Allura XPER FD10/FD20		12NC: 4522 901 00140		Version: 00 Doc.Date: 2006-09-01	
Author: W.Bleyen		Pag.: 3 of 3	Group No: 579	Section: -	Status: Approved