


Product name:	AC Supply for HeartStart MRx Defibrillator	
Identification code	M3539A	
Total weight (in Kg)	Up to .8 kgs	
Producer /Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA
	Electronic info:	www.philips.pms.recycleinfo.com

Recycle info	Substances:	Location
Special attention	Item	Location
	After removing 4 screws, one of which is under a label, pry case along seam. Open to locate pc board inside.	

Batteries	Type:	Location
	None	
To be removed		

Hazardous	None	
		
To be removed		

Fluids / Gases	Item	Location
	None	