


Product name	Invivo Sense Neuro Vascular coils
Identification code(s)	4522 132 31490 8-SNVC 1.5T 4522 132 31710 8-SNVC 3.0T 4522 132 31080 16-SNVC 1.5T 4522 132 31720 16-SNVC 3.0T 9896 030 13370 4-NV 3.0T
Total weight [kg]	SNVC → 17 Kg 4-NV → 12 Kg
Manufacturer	Invivo, corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_NeuroVascularCoils.doc

Recycling information

Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure


Hazardous substances	Items	Location
	Several printed circuit boards inside the coils	Figures next page
	Beryllium copper contact pins on connectors	Figures next page

Name: Recycling passport Invivo SNVC				12NC: 4522 132 31490		Version: 01 Doc.Date: 2007-05-11	
Author: Paul Ebben		Page 1 of 3		Group No: 5790	Section: -	Status: Authorized	

Fluids and gases	Items	Location
	Plastic bottle with copper sulfate (4-NV coil only)	Figure next page
	0.008moles/liter CuSO ₄ and 0.034moles/liter	
	Bottle is delivered with the coil separately	

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



Name: Recycling passport Invivo SNVC	12NC: 4522 132 31490	Version: 01 Doc.Date: 2007-05-11
Author: Paul Ebben	Page 2 of 3	Group No: 5790 Section: -
Status: Authorized		

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

Written By Paul Ebben	Reviewed By (optional) Geert-Jan Platel
Approved By Paul Ebben	Approved By (optional)

End of document.

Name: Recycling passport Invivo SNVC	12NC: 4522 132 31490	Version: 01 Doc.Date: 2007-05-11
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: - Status: Authorized