


Product name	System Filter Box
Identification code(s)	4522 132 5600 x 4522 132 5730 x 4522 132 2465 x 4522 131 3421 x 4522 131 6089 x
Total weight [kg]	110 kg
Manufacturer	INNO/Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	System_SystemFilterBox_452213224650_579_v03.doc

For picture of component : see figure 1


Recycling information

Special Attention 	Items	Location
	None	
Batteries 	Items, type	Location
Hazardous substances 	Items	Location
	300V power supply	See figure 2
	Astec power supply	See figure 2
	Combifilter	See figure 1

Name: Recycling Passport SystemFilterBox			12NC: 4522 132 24650			Version: 03 Doc.Date: 2007-08-01		
Author: J. de Haas			Pag.: 1 of 3	Group No: 5790	Section: -		Status: Authorized	

Locations as mentioned in the Passport (Information in pictures)

Figure 1




Name: Recycling Passport SystemFilterBox				12NC: 4522 132 24650		Version: 03 Doc.Date: 2007-08-01	
Author: J. de Haas		Pag.: 2 of 3		Group No: 5790	Section: -	Status: Authorized	

© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007

Template version 2.

Figure 2



Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2006-07-26
2	Updated version	2007-05-20

Recycling Passport Approval Signatures

Written By	Reviewed By (optional)
J.de Haas	
Approved By	Approved By
GJ.Plattel	

End of document.

Name:	Recycling Passport SystemFilterBox			12NC:	4522 132 24650			Version: 03
Author: J. de Haas			Pag.: 3 of 3	Group No: 5790	Section: -		Status: Authorized	Doc.Date: 2007-08-01