

|                        |                                                                |
|------------------------|----------------------------------------------------------------|
| <b>Product name</b>    | <b>NT-DAC 1.5T-C-XW6200 R1</b>                                 |
| Identification code(s) | 4522 132 5076x                                                 |
| Total weight [kg]      | 137 kg excluding CDAS ( $\pm 27$ kg) and XW6200 ( $\pm 17$ kg) |
| <b>Manufacturer</b>    | <b>Philips Medical Systems Nederland B.V.</b>                  |
| Address                | Veenpluis 4-6 Best                                             |
| Zip code               | 5680DA                                                         |
| Country                | The Netherlands                                                |
| Electronic info        | System_NTDAC1.5TforCDAS_452213250760_579_v01.doc               |

## Recycling information

| Special Attention                                                                 | Items                                            | Location |
|-----------------------------------------------------------------------------------|--------------------------------------------------|----------|
|  | <b>NOT part of this passport:</b>                |          |
|                                                                                   | 1/ CDAS (45221320893x)                           | Figure 2 |
|                                                                                   | 2/ XW6200(45221324545x)                          | Figure 3 |
|                                                                                   | 3/ Stand. Cab. 33HE + 2HE Lade (4522 131 7443 x) | Figure 1 |
|                                                                                   |                                                  |          |

| Batteries                                                                           | Items, type    | Location |
|-------------------------------------------------------------------------------------|----------------|----------|
|  | Not applicable |          |
|                                                                                     |                |          |
|                                                                                     |                |          |
|                                                                                     |                |          |

| Hazardous substances                                                                | Items                                           | Location |
|-------------------------------------------------------------------------------------|-------------------------------------------------|----------|
|  | Cable bracket front (Zn yellow passivated 5µm)  | Figure 2 |
|                                                                                     | (quick acting) Screw (Zn yellow passivated)     | Figure 3 |
|                                                                                     | Cable bracket back (Zn yellow passivated 5µm)   | Figure 5 |
|                                                                                     | Cable bracket (Zn yellow passivated 5µm)        | Figure 6 |
|                                                                                     | Cable relief bracket (Zn yellow passivated 5µm) | Figure 7 |
|                                                                                     |                                                 |          |

| Fluids and gases                                                                    | Items          | Location |
|-------------------------------------------------------------------------------------|----------------|----------|
|  | Not applicable |          |
|                                                                                     |                |          |
|                                                                                     |                |          |
|                                                                                     |                |          |

|                                                            |                                |                                                              |
|------------------------------------------------------------|--------------------------------|--------------------------------------------------------------|
| Name:<br><b>Recycling passport NT-DAC 1.5T-C-XW6200 R1</b> | 12NC:<br><b>4522 132 50760</b> | Version: <b>01</b><br>Doc.Date: <b>2007-05-11</b>            |
| Author: <b>G. Geboers</b>                                  | Page <b>1</b> of <b>3</b>      | Group No: <b>5790</b> Section: <b>-</b> Status: <b>Auth.</b> |

## Locations as mentioned in the Passport (Information in pictures)



Cdas



Figure 1 Frontview  
NT-DAC cabinet



Stand. Cabinet

Figure 2



Figure 3



XW 6200




Figure 5



Figure 6

Figure 7



Figure 4 Backview  
NT-DAC cabinet



|                                                            |  |  |                                |                       |                   |                                                   |                      |  |
|------------------------------------------------------------|--|--|--------------------------------|-----------------------|-------------------|---------------------------------------------------|----------------------|--|
| Name:<br><b>Recycling passport NT-DAC 1.5T-C-XW6200 R1</b> |  |  | 12NC:<br><b>4522 132 50760</b> |                       |                   | Version: <b>01</b><br>Doc.Date: <b>2007-05-11</b> |                      |  |
| Author: <b>G. Geboers</b>                                  |  |  | Page <b>2</b> of <b>3</b>      | Group No: <b>5790</b> | Section: <b>-</b> |                                                   | Status: <b>Auth.</b> |  |

## Recycling Passport Revision History

| Revision | Description     | Date       |
|----------|-----------------|------------|
| 1        | Initial version | 2006-08-18 |

## Recycling Passport Approval Signatures

|                                             |                               |
|---------------------------------------------|-------------------------------|
| <b>Written By</b><br><br>G. Geboers         | <b>Reviewed By (optional)</b> |
| <b>Approved By</b><br><br>Geert-Jan Plattèl | <b>Approved By</b>            |

End of document.

|                                                            |                                |                                                              |
|------------------------------------------------------------|--------------------------------|--------------------------------------------------------------|
| Name:<br><b>Recycling passport NT-DAC 1.5T-C-XW6200 R1</b> | 12NC:<br><b>4522 132 50760</b> | Version: <b>01</b><br>Doc.Date: <b>2007-05-11</b>            |
| Author: <b>G. Geboers</b>                                  | Page <b>3</b> of <b>3</b>      | Group No: <b>5790</b> Section: <b>-</b> Status: <b>Auth.</b> |