Philips Medical Systems/ Product Recycling Passport

Product name: 6-S		6-Slot N	Aodule Rack	
Identification code(s)		M1276A		
Total weight (in Kg)		1.1		
Producer/	Name company:		Philips Medical Systems	
Manufacturer	Address:		3000 Minuteman Road, Andover, MA	
	Zip code:		01810	
	Country:		USA	

Recycle Info	Items:	Location
Special attention	Separate the front assembly from the rear housing assembly.	Figure 1
Fluids / Gases	Items: Not Applicable	Location
Batteries	Type: Not Applicable	Location
To be Removed		
Hazardous To be Removed	Substances: Printed circuit board, lead solder (Pb)	Location Figure 2



Title: M1276A Satellite Recycling Passport Num

Philips Medical Systems/ Product Recycling Passport






Figure 2

Title: M1276A Satellite 6-Slot Module Rack Recycling Passport Recycling Passport Number: A-M1276-90021 Rev: A