

Product name:	<b>7T Achieva MR Scanner RF Cabinet</b>	
Identification code	9896 030 5059*	
Total weight (in Kg)	Approximately 485 Kg	
Producer /Manufacturer	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/mri/">http://www.medical.philips.com/main/company/sustainability/recycling/mri/</a>

Recycle info	Substances:	Location
<b>Hazardous</b>  <b>To be removed</b>	Lead (Pb) (Solder - PC boards)	Figure 1
 <b>Batteries</b> <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	Not applicable.	
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
	Storage capacitors	Figure 1,2
 <b>Fluids / Gases</b>	<b>Item</b>	<b>Location</b>
	Not applicable.	

\* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport ( information)



Pictures




Figure 1. 7T RF cabinet, front view.



Figure 2. 7T RF Cabinet, rear view.