**Recycling Passport** 

Product name:	DigitalDiagnost THS			
Identification code	9890-010-8	4301		
Order No.				
Device Type	B2B (Business to business)			
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH		
	Address:	Roentgenstrasse 24		
	Zip code:	22335 Hamburg		
	Country:	Germany		
	Web info:	See Philips Intranet		
Total weight (in Kg)	•		212 kg	
thereof	•			

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	98,244	KG
	iron, high alloy (>5%)	34,198	KG
steel, iron		132,442	KG
nonferrous metals and alloys	aluminum, -alloy	38,655	KG
	copper, -alloy	0,047	KG
nonferrous metals and alloys		38,702	KG
plastics / organic substances	thermoset	17,163	KG
	thermoplastic	7,376	KG
	hydraulic oil	2	KG
	elastomer	0,018	KG
plastics / organic substances		26,557	KG
standard parts	other electronic powered devic	5,224	KG
	motors, pumps	4	KG
	mounting parts, attaching part	3,423	KG
	cables	1,388	KG
	printed circuit boards	0,402	KG
standard parts		14,437	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	confirmed by supplier		
prohibited materials			
		212,138	KG

Printed copies are uncontrolled documents


Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements.  For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous	Substances:	Location
Substances		200000
To be removed		
Batteries	Type:	Location
To		2300000
_		
be removed		
Special attention	Item	Location
	Before dismantling the hydraulic oil, locked the body of the column with the two red bolt.	Fig. 2
Eluida / Corre		T
Fluids / Gases	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1
	HV-transformer oil contains no PCBs	
	•	

SAP-Doc.:9890-010-84301   Doctype: 129-M03   SAP-Vers.: 00   Date : 22.10	22.10.2010	Page2/3
---	------------	---------

## recycling i asspo

 $Locations \ as \ mentioned \ in \ the \ Passport \ \ (\ Pictures \ information)$ 





Printed copies are uncontrolled documents

SAP-Doc.:9890-010-84301 | Doc.-type: 129-M03 | SAP-Vers.: 00 | Date : 22.10.2010 | Page3/3

## **Recycling Passport**

Printed copies are uncontrolled documents

SAP-Doc.:9890-010-84301 | Doc.-type: 129-M03 | SAP-Vers.: 00 | Date : 22.10.2010 | Page4/3