

<b>Product name:</b>	<b>7T Achieva MRI Nova TR Coil + Nova 16Ch Receive coil</b>	
<b>Identification code</b>	4535 674 2804* TR Coil 4535 674 2805* 16 Channel Receive Coil	
<b>Total weight (in Kg)</b>	Approximately 13.0 Kg	
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/mri/">http://www.medical.philips.com/main/company/sustainability/recycling/mri/</a>

Recycle info	<b>Substances:</b>	Location
 <b>Hazardous</b>  <b>To be removed</b>	Lead(Pb) (Solder – PC boards)	Figure 1, 2
 <b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	Not applicable.	
 <b>Special attention</b>	<b>Item</b>	<b>Location</b>
	Not applicable	
 <b>Solids / Gases</b>	<b>Item</b>	<b>Location</b>
	Not applicable	

\* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)



Figure 1. 7T Achieva Nova TR Coil, and Nova 16 Channel receive coil, front view



Figure 2. 7T Achieva Nova TR Coil, and Nova 16 Channel receive coil, rear view.



Figure 3. Nova TR and 16 Channel coils in nested configuration. Rear and front views.