

Philips Medical Systems / Product Recycling Passport

Product name:	Anesthetic Gas Module	
Identification code	M1026B	
Total weight (in Kg)	6.40 kg	
Producer /Manufacturer	Name company:	Philips Medizin Systeme Boeblingen
	Address:	Hewlett-Packard Strasse 2
	Zip code:	Boeblingen 71034
	Country:	Germany

Recycle info	Substances:	Location
Hazardous  To be removed		
Batteries  To be removed	Type:	Location
	Not applicable.	
Special attention 	Item	Location
Biohazard 	Item	Location
	Bio-contaminated tubings + filters	Figure 1-3

Material content (global)	Weight in kg	Materials content (continue)	Weight in kg
Total PCB	0.63		
Plastics	0.53		
Steel	4.91		
Aluminum	0.33		

Locations as mentioned in the Passport (Pictures information)




Figure 1




Figure 2




Figure 3