Recycling Passport


Product name:	Optimus 30/40		
Identification code	9890-010-82591		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
Total weight (in Kg)			150 Kg
thereof			

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	61,63	KG
steel, iron		61,63	KG
nonferrous metals and alloys	aluminum, -alloy	2,142	KG
nonferrous metals and alloys		2,142	KG
glass / ceramics	ferrite	0,5	KG
glass / ceramics		0,5	KG
plastics / organic substances	oil	4	KG
	thermoplastic	1,2	KG
	thermoset	0,5	KG
plastics / organic substances		5,7	KG
standard parts	power supplies / transformers	27,881	KG
	cables	24,711	KG
	other electronic powered devic	13,9	KG
	printed circuit boards	10	KG
	mounting parts, attaching part	3,295	KG
standard parts		79,787	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	confirmed by supplier		
prohibited materials			
		149,759	KG

Printed copies are uncontrolled documents

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous	Substances:	Location
To be removed		
Batteries	Type:	Location
To be removed		Location
Special attention	Item	Location
		Zocuron
Fluids / Gases	Item	Location
	Oil tank / Remove oil before dismantling HV-transformer oil contains no PCBs	Fig. 1

Locations as mentioned in the Passport (Pictures information)



Printed copies are uncontrolled documents

SAP-Doc.:9890-010-82591 | Doc.-type: 129-M03 | SAP-Vers.: 00 | Date. 22.10.2010 | Page3/3