
Product name: il		iU22 M	iU22 Medical Diagnostic Ultrasound Equipment	
Identification code(s)		8500-0064-02		
Total weight (in Kg)		160 Kg		
Dece dece ere/	Name company: Address:		Philips Ultrasound	
Producer/			Filinps Oldasound	
Manufacturer			22100 Bothell Everett Highway	
	Zip code:		98041-3003	
	Country:		USA	

Recycle Info	Items:	Location
Special attention	Gas Spring under pressure. See Figure 1 [1]	
Fluids / Gases	Items:	Location
	Oil and Nitrogen Gas: (2X) in Gas Spring Tubes in Control Panel Articulation Assy Do not drain oil with device under pressure.	Figure 1 [1]
	Oil and Nitrogen Gas: in Gas Spring Tube in Arm Articulation Assy Do not drain oil with device under pressure.	Figure 1 [2]
Batteries	Туре:	Location
	Li Battery: (2X) on CPU PCB	Figure 1 [3]
	Battery Lead Acid (Pb): in Power Supply Assy	Figure 1 [4]
To be Removed		
Hazardous	Substances:	Location
To be Removed	Lead (Pb) - contains lead in PCB's, solder, wiring, and cables	N/A
	Mercury (Hg) – Video monitor lamps	Figure 2 [5]
	Mercury (Hg) – Touch screen lamps	Figure 2 [6]
	Note: Transducers (Probes) contain lead (Pb) on Flex circuits / PCB's in the connector and probe areas	N/A


Philips Medical Systems/ Product Recycling Passport





Title: iU22 Recycling Passport Recycling Passport Number: 8500-0064-00050 Rev: B





Recycling Passport Revision History

Re	ev	Description	Date
A	4	First Passport	July 31, 2006
В	3	Revision for Product Weight Correction	September 14, 2009

Title: iU22 Recycling Passport Recycling Passport Number: 8500-0064-00050 Rev: B