




PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 2
-------------------------------	-----------------------------------	-------------

Product name:	PFEI-C-8CH-30T	
Identification code(s)	4522-133-0438x	
Total weight (in Kg)	20 Kg	
Producer /Manufacturer	Name company:	Philips Medical Systems Best
	Address:	Veenpluis 4-6
	Zip code:	5680DA Best
	Country:	The Netherlands
	Electronic info:	System_PFEI-C-8ch-30T_452213304380_579_v00.doc

Recycle info	Substances:	Location
<i>Special attention</i> 	Items:	
	No special attention needed.	
 Batteries <i>To be removed</i>	Type:	Location
	No batteries present.	
<i>Hazardous</i>  <i>To be removed</i>	Examples	Location
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Locations as mentioned in the Passport (Pictures information)




Figure 1a, PFEI-C-30T right-side




Figure 1b, PFEI-C-30T left-side




Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly