

Product name:	3rd Order Static Shim Tube/Coils		
Identification code	4535 670 5747*		
Total weight (in Kg)	Approximately 200 Kg		
Producer /Manufacturer	Name company:	Philips Medical Systems (Cleveland), Inc	
	Address:	595 Miner Rd., Cleveland, Ohio	
	Zip code:	44143	
	Country:	USA	
	Electronic info:	http://www.medical.philips.com/main/company/sustainability/recycling/mri/	

Recycle info	Substances	Location
Hazardous  To be removed	None.	
 Batteries To be removed	Type	Location
	None.	
Special attention 	Item	Location
	None.	
 Fluids / Gases	Item	Location
	None.	

* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)



Figure 11 7T 3rd order shim tube, built with lead-free (silver) solder. Illustrated as removed from magnet.