

PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 3
-------------------------------	-----------------------------------	-------------

Product name:	Brilliance CT Reconstruction Cabinet	
Identification code(s)	4535 674 12261, 4535 674 12271, 4535 674 12031, 4535 674 65881, 4535 674 12451, 4535 674 12461, 4535 674 12471	
Total weight (in Kg)	Maximum possible weight of 227 Kg. Configuration dependent.	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	595 Miner Road Highland Heights, Ohio
	Zip code:	44143
	Country:	United States of America

Printed copies are uncontrolled

Recycle Info	Items:	Location
Special attention 		
Fluids / Gases 		
Batteries  To be Removed	Type: NiCd batteries (typical location)	Location Figure 1 (2)
Hazardous  To be Removed	Substances: Printed Circuit Boards (typical location)	Location Figure 1 (1)

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




Figure 1: CIRS Rack

PHILIPS HEALTHCARE	Product Recycling Passport	Page 3 of 3
-------------------------------	-----------------------------------	-------------

Recycling Passport Revision History

Rev	Description	Date
C	CO 2015923 Release of online WEEE passport	2 May 2007
D	CO 2019848 – Changed Identification label to Product Part Number	25 February 2009
E	CO 2020759 - Removed superfluous footer, added additional part number to identification code	9 June 2009

Printed copies are uncontrolled