

<b>PHILIPS HEALTHCARE</b>	<b>Product Recycling Passport</b>	Page 1 of 3
-------------------------------	-----------------------------------	-------------


<b>Product name:</b>	<b>Brilliance Mobile Network Device (MND)</b>	
<b>Identification code(s)</b>	<b>455011002301 and 455011002311</b>	
<b>Total weight (in Kg)</b>	<b>15 Kg (455011002301), 50 Kg (455011002311)</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems
	Address:	595 Miner Road Highland Heights, Ohio
	Zip code:	44143
	Country:	United States of America

Printed copies are uncontrolled

Recycle Info	Items:	Location
<b>Special attention</b> 		
<b>Fluids / Gases</b> 		
<b>Batteries</b>  <b>To be Removed</b>	<b>Type:</b>	<b>Location</b>
	<b>Li batteries (typical location)</b>	Figure 1 (1)
<b>Hazardous</b>  <b>To be Removed</b>	<b>Substances:</b>	<b>Location</b>
	<b>Printed Circuit Boards (typical location)</b>	Figure 1 (2)

**Figure 1: Show locations of materials mentioned on the previous sheet. (Picture information)**


**Brilliance Mobile Network Device (MND) Part Number 455011002301**



**Brilliance Mobile Network Device (MND) Part Number 455011002311**

Tower option

Rack mount option



THE: recycling passport of Brilliance Mobile Network Device (MND)  
Recycling Passport Number: 455019702431 Rev: A

Printed copies are uncontrolled

### Recycling Passport Revision History

Rev	Description	Date
A	CO H3002318 Release	March 26, 2009

Printed copies are uncontrolled