


Product name	Philips 3.0T 32-Channel Digital Head Coil
Identification code(s)	9896-032-06931
Total weight [kg]	10.43 kg
Manufacturer	Invivo
Address	3545 SW 47 th Ave., Gainesville, FL
Zip code	32608
Country	USA
Electronic info	208161.doc

Recycling information

Special Attention	Items	Location
	None	N/A

Batteries	Items, type	Location
	None	N/A

Hazardous substances	Items	Location
	Circuit boards within the coil with lead solder	Figure 1
	Beryllium copper on contact pins on connectors	
	and three springs in anterior latch assemblies	

Fluids and gases	Items	Location
	None	N/A

Name:		12NC:		Doc.Date: 5/23/2012	
Recycling Passport Philips 3.0T 32-CH Digital Head Coil		9896-302-06931			
	Pag.: 1 of 2	Group No:	Section: -	Status: Rev 1	

Locations as mentioned in the Passport (Information in pictures)



Figure 1

Name: Recycling Passport Philips 3.0T 32-CH Digital Head Coil				12NC: 9896-302-06931		Doc.Date: 5/23/2012	
			Pag.: 2 of 2	Group No:	Section: -	Status: Rev 1	
© PHILIPS HEALTHCARE NEDERLAND BV, BEST 2008				Ref. XJV032011.07 Rev.01			