


Recycling Passport

Product name:	OPTIMUS CXA 80 1T 480V		
Identification code	9890-000-70601; 9890-000-70605		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.healthcare.philips.com/main/about/sustainability/recycling/xraygenerators.wpd	
Total weight (in Kg)			268
thereof			

steel, iron	iron, low alloy (<5%)	113,568	KG
	iron, high alloy (>5%)	0,162	KG
steel, iron		113,73	KG
nonferrous metals and alloys	copper, -alloy	3,96	KG
	aluminum, -alloy	1,103	KG
	zinc, -alloy	0,883	KG
nonferrous metals and alloys		5,946	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	53,1	KG
	thermoplastic	5,663	KG
	thermoset	4,558	KG
	elastomer	0,247	KG
plastics / organic substances		63,568	KG
standard parts	other electronic powered devic	61,769	KG
	printed circuit boards	10,52	KG
	cables	9,356	KG
	mounting parts, attaching part	2,846	KG
	motors, pumps	0,276	KG
	fans	0,064	KG
standard parts		84,831	KG
		268,363	KG

Additional weight for Systemcomponents see page three.

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances 	Substances:	Location
To be removed		
Batteries 	Type:	Location
	Lithium chrome cell, 3V	Figure 1 (1)
	Lithium chrome cell, 3V	Figure 1 (2)
To be removed		
Special attention 	Item	Location
	Pb is present in the soldering process of PCBs	Figure 1, 2
Fluids / Gases 	Item	Location
	Transformer Oil, Shell DIALA	Figure 1 (3)
	HV-transformer oil contains no PCBs	

Recycling Passport




Figure 1

4. Room in/out Code Nr. 4512-201-0417. overall weight 1,9 kg. Materials see Recycling Passport.

4.1 Cablesset cabinet Code Nr. 4512-201-0567. overall weight 0,65 kg. Materials see Recycling passport.

4.2 System Power distribution Unit Code Nr. 4512-201-0476. overall weight 4,8 kg. Materials see Recycling Passport.