





Recycling Passport

Product Name:		CT Scanner System
Product Model		MX 16-slice
Identification Code		989605656511
Total weight (in Kg)		2210
Producer /Manufacturer	Name	Philips and Neusoft Medical Systems Co., Ltd
	Address:	Neusoft Park, Hun Nan Industrial Area, Shenyang 110179, China

Recycle info	Substances:	Location
Hazardous  To be removed	Examples	
	Lead (Pb)	Figure 2
Batteries  To be removed	Type:	Location
	Lithium Coin Battery	Figure 5
		Figure 10
Special attention 	Item	Location
	Air-spring	Figure 1
		Figure 7
Fluids / Gases 	Item	Location
	High pressure oil tank	Figure 2

MX 16-slice:

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Lead (Pb)	25.151	Printed Circuits Boards	6.874
Iron (Fe)	1259.57	Tungsten (W)	2.183
Aluminum (Al)	563.69		
Copper (Cu)	15.828	All other material types	336.704

Locations as mentioned in the Passport (Pictures information)




Figure 1 CT system front view




Figure 2 Gantry front views (without cover)




Figure 3 Gantry front views (without cover)




Figure 4 Gantry rear view (without cover)




Figure 5 Gantry Left side (without cover)




Figure 6 Gantry right side (without cover)




Figure 7 Couch System (without cover)




Figure 8 Couch (inside view)




Figure 9 Computer case front view




Figure 10 Computer case front view (door open)




Figure 11 Computer case rear view