## Product (MX40) Recycling Passport

Notice UNCONTROLLED IN PRINTED FORM UNLESS STAMPED IN RED

**Document Ownership**  The document owner is listed below:

Function	Title
R&D	Project Leader

The current document revision is described below:

Revision Summary

Rev.	Change Order #	Change Summary
Α	Enter EC number	As Issued

## Product (MX40) Recycling Passport

Product nar	ne:	Enter name of product as it appears on the product. MX40 Patient Worn Monitor		
Reference Ider	ntifier:	<i>This should be the model number as it appears on the product.</i> 865350, 865351, 865352		
Total weight (in Kg):0.403		0.403 1	3 Kg	
Producer /Manufacturer			Philips Healthcare 3000 Minuteman Road	
	Zip code:		01810	
	Country:		USA	

Recycle info	Substances:	Location
Special attention	Indicate any part of the product that could pose a hazard to the recycler and treatment facility during disassembly such as capacitors with residual energy or spring loaded valves and show location	
	TFT-Display	Figure 1 & Figure 4
Batteries	Type:	Location
2	Indicate all batteries in the product via battery chemistry and show location	
	Lithium Ion battery (optional) Lithium Ion battery (coin cell), not removable	Figure 2 Figure 5
To be removed		
Hazardous	Item	Location
	Indicate location of printed circuit boards and mercury in back lights as well as any cadmium, flame retardants (PBB and PBB) and hexavalent chromium in the product and show location	
3K	See table below	
To be removed		

## Additional Information as needed

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Battery	~0.1		
Plastics etc	~0.2		
Printed Circuits Boards	~0.1		

## Product (MX40) Recycling Passport

Locations as mentioned in the Passport (Pictures information)

