


Recycling Passport

Product name:	Velara R/F 2T	Velara R/F 2T	Velara R/F Rel.2
Identification code	9890-000-62062	9890-000-62063	9890-000-70151
Product name:	Velara R/F 2T Rel.3	Velara R/F 2T Rel.3	
Identification code	9890-000-70291	9890-000-70294	
Product name:	Velara R/F 2T R.3.11.1	Velara R/F 2T R.3.12.1	
Identification code	9890-000-70295	9890-000-70296, 9890-000-70297	
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.medical.philips.com/main/company/sustainability/recycling/xray_generators/	
Total weight (in Kg)			285
thereof			

steel, iron	iron, low alloy (<5%)	130,785	KG
	iron, high alloy (>5%)	2,334	KG
steel, iron		133,119	KG
nonferrous metals and alloys	aluminum, -alloy	12,274	KG
	copper, -alloy	5,005	KG
	zinc, -alloy	1,576	KG
nonferrous metals and alloys		18,855	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	48,06	KG
	thermoplastic	5,952	KG
	thermoset	4,245	KG
	elastomer	0,232	KG
plastics / organic substances		58,489	KG
standard parts	other electronic powered devices	45,92	KG
	printed circuit boards	11,584	KG
	cables	6,733	KG
	power supplies / transformers	6,184	KG
	mounting parts, attaching part	2,03	KG
	fans	1,804	KG
	motors, pumps	0,52	KG
standard parts		74,775	KG
		285,526	KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances 	Substances:	Location
To be removed		
Batteries 	Type:	Location
	Lithium chrome cell, 3V	Figure 1 (1)
	Lithium chrome cell, 3V	Figure 1 (2)
To be removed		
Special attention 	Item	Location
	Pb is present in the soldering process of PCBs	Figure 1
Fluids / Gases 	Item	Location
	Transformer Oil, Shell DIALA	Figure 1 (3)
	HV-transformer oil contains no PCBs	

Recycling Passport




Figure 1