





Product name:	PFEI-C-10T		
Identification code(s)	4522-132-2023		
Total weight (in Kg)	20 Kg		
Producer /Manufacturer	Name company:	Philips Medical Systems Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	PFEI-C-10T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed	Examples	Location
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: Recycling Passport PFEI-C-10T		12NC: 4522 132 20230		Version: 04 Doc.Date: 2007-08-08	
Author:		Pag.: 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Pictures information)




Figure 1a, PFEI-C-10T right-side




Figure 1b, PFEI-C-10T left-side

Name: Recycling Passport PFEI-C-10T		12NC: 4522 132 20230		Version: 04 Doc.Date: 2007-08-08	
Author:		Pag.: 2 of 3	Group No: 5790	Section: -	Status: Authorized
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007				Template version 2.	




Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly

Name: Recycling Passport PFEI-C-10T		12NC: 4522 132 20230		Version: 04 Doc.Date: 2007-08-08	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: Authorized
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007				Template version 2.	