





Product name	NT-DAC
Identification code(s)	98960302243x / 98960302227x / 98960302529x / 98960302528x
Total weight [kg]	± 137 kg excluding CDAS (±27kg)/ Z420 (+/-14Kg) / Z400 (±17kg) and Switch (±2kg)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

Recycling information

Special Attention 	Items: The following items (if applicable) have their own passport and are not part of this passport.	Location
	1/ STAND.CAB. 33E 45221317443x	Figure 1
	2/ CDAS II 45980001447x	Figure 1
	3/ Z420 45980032730x / 45980032733x Z400 45221331444x	Figure 1
	4/ SWITCH (HP Procurve) 45980001377x Edgecore 45980020324x	Figure 1


Batteries 	Items, type	Location
	Not applicable	

Hazardous substances 	Items	Location
	Cable bracket front (Zn yellow passivated 5µm)	Figure 1
	Cable bracket back (Zn yellow passivated 5µm)	Figure 2
	Cable bracket (Zn yellow passivated 5µm)	Figure 2
	Cable relief bracket (Zn yellow passivated 5µm)	Figure 2


Fluids and gases 	Items	Location
	Not applicable	

Name: WEEE RECYCLING PASSPORT NT-DAC	12NC: DMR164513	Version: 01 Doc.Date: 2013-04-18
Author:	Pag.: 1 of 2	Group No: 5790 Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)



Figuur 1 Frontview
NT-DAC



Figuur 2 Rearview
NT-DAC

Name: WEEE RECYCLING PASSPORT NT-DAC		12NC: DMR164513	Version: 01 Doc.Date: 2013-04-18
Author:	Pag.: 2 of 2	Group No: 5790	Section: -
Status: Authorized			Template version 2.