




Product name	CDAS (PFS-316) (1.5T), CDAS (PFS-316) (3T)		
Identification code(s)	4522 132 0893x	4522 132 5996x	4598 000 1447x
	4522 132 0894x	4522 133 0471x	
Total weight [kg]	33.6 kg (includes packaging) / 26 kg (without packaging)		
Manufacturer	RadiSys Corporation		
Address	5445 NE Dawson Creek Drive, Hillsboro, OR		
Zip code	97124		
Country	USA		
Electronic info			

Recycling information

Special Attention	Items	Location
	Power Supply, purchased as sealed unit, contents uncertain.	Figure 1


Batteries	Items, type	Location
	CPU board Coin Battery, Lithium	Figure 1

Hazardous substances	Items	Location
	CPU PWB (Remark: PWB = Printed Circuit Board)	Figure 1
	Other PWB's present in other slots (incl. emigaskets)	Not shown
	Fan Controller PWB	Figure 1
	Backplane, Main PWB	Figure 1
	Backplane, R(ear)T(ransission)M(odule) PWB	Figure 1
	EMI gaskets (contact springs) on front- /rearplates / fan tray cover and dustfilter cover contain BeCu	Figure 1
	Power Supply <i>may</i> contain: PWB, Electrolyte Capacitors	Figure 1

Name: Recycling passport CDAS ASSY II	12NC: DMR108299	Version: 00 Doc.Date: 2009-11-27
Author:	Page 1 of 2	Group No: 5790 Section: -
		Status: Approve

Locations as mentioned in the Passport (Information in pictures)

Figure 1: Front of CDAS Chassis Assembly



Name: Recycling passport CDAS ASSY II	12NC: DMR108299	Version: 00 Doc.Date: 2009-11-27
Author:	Page 2 of 2	Group No: 5790 Section: -
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2010		Status: Approve