




Product name	BODY COIL WA 3.0T
Identification code(s)	9896 030 2173.0
Total weight [kg]	60 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

Recycling information

Special Attention	Items	Location
	NA	Figure ..


Batteries	Items, type	Location
	NA	Figure ..

Hazardous substances	Items:	Location, see below
	4522 117 9840.X RF SCREEN PART A (2X)	1
	4522 117 9841.X RF SCREEN PART B (2X)	2
	4522 132 8043.X MATCHING GROUND FLEX (2X)	3
	4522 132 7896.X CONNECT BRACKET PCB (1X)	4
	4522 117 9825.X DETUNE BOARD 3.0T (16X)	5
	4522 132 8044.X BC SOLDER FLEX RF SCREEN (32X)	6
	4522 132 8045.X BC SOLDER FLEX ANTENNA (4X)	7
	4522 117 9829.X MATCHING BOARD 3.0T (2X)	8
	4522 133 1045.X CABLE MZ-X1- MVM-X1 (MIC) (1X)	9

Name: Recycling Passport BODY COIL WA 3.0T	12NC: 9896 030 2173.0	Version: 00 Doc.Date: 2010-11-15
Author:	Pag.: 1 of 2	Status: Authorized

Fluids and gases	Items	Location
	NA	Figure ..

Locations as mentioned in the Passport (Information in pictures)



Name: Recycling Passport BODY COIL WA 3.0T		12NC: 9896 030 2173.0		Version: 00 Doc.Date: 2010-11-15	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: Authorized
© PHILIPS HEALTHCARE NEDERLAND BV, BEST 2008			Ref. XJV032011.07 Rev.01		