



Product name	Invivo Wrist coils	
Identification code(s)	4522 150 40360 4-Wrist 1.0T 9896 030 13420 4-Wrist 1.5T 9896 030 13390 4-Wrist 3.0T	4522 132 63910 8-Wrist 1.5T 4522 132 63920 8-Wrist 3.0T
Total weight [kg]	4-ch Wrist 6 kg 8-ch Wrist 7 Kg	
Manufacturer	Invivo, corporation. USA	
Address	3545 SW 47th Avenue	
Zip code	Gainesville, FL 32608	
Country	USA	
Electronic info	RFCoils_Invivo_WristCoils V2.doc	

Recycling information

Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure

Hazardous substances	Items	Location
	Beryllium copper contact pins on connectors	Figures next page
	Several printed circuit boards inside the coils	Figures next page

Fluids and gases	Items	Location
	4-ch coils; Plastic bottle with copper sulfate 0.008moles/liter CuSO4 and 0.034moles/liter NaCl	Figures next page
	8-ch coils; Plastic bottle with Nickel Chloride 31.5mM/L of NiCL2-6H2O	Figures next page

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:




Figure 1 Wrist 4-ch




Figure 2, Wrist 4-ch





Figure 3 Wrist 8-ch