





<b>Achieva</b>	<b>Achieva 3.0T TX Covers</b>		
Identification code(s)	SET COVERS 3.0T TX :		
	MAGNET FRONT COVERS CYL-A-TX	9896 030 2125x	
	MAGNET REAR COVERS CYL-A-TX	9896 030 2131x	
Total weight [kg]	± 150 kg		
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>		
Address	Veenpluis 4-6 Best		
Zip code	5680DA		
Country	The Netherlands		
Electronic info			

### Recycling information

Special attention	Items	Location
	Products are made of composite materials consisting of glass, fibers and resins.	
	Microphone printed circuit board in front cover	Figure 1
	4x light sources (LED printed circuit board)	Figure 1


Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Not applicable	

Fluids and gases	Items	Location
	Not applicable	

Name:	<b>Recycling passport Covers 3.0T TX</b>	12NC:	<b>9896 030 21250</b>	Version: <b>00</b>
Author:		Group No: 5790	Section: -	Doc.Date: <b>2009-07-17</b>
	Pag.: <b>1 of 2</b>			Status: <b>Authorized</b>

**Overview of covers:**



**figure 1**

Name:	<b>Recycling passport Covers 3.0T TX</b>	12NC:	<b>9896 030 21250</b>	Version: <b>00</b>	Doc.Date: <b>2009-07-17</b>
Author:		Pag.: <b>2 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>