



**UNIKLINIK
KÖLN**

Institut und Poliklinik für
Radiologische Diagnostik




Clinical Benefit of Spectral CT

Medical Imaging Convention London ExCel 2018 |

David Maintz



Computed Tomography (CT)



1972: EMI CAT Scanner







Spectral CT today






Background CT




Background Spectral CT



- Attenuation μ behaves differently for different materials at different x-ray energies → CT numbers depend on energy
- Only water is “forced” to 0 HU


$$\text{CT}_{\text{pixel}} = \frac{\mu_{\text{pixel}} - \mu_{\text{water}}}{\mu_{\text{water}}} * 1000 \text{ HU}$$

Interaction of photons with matter



Background Spectral CT

Interaction of photons with matter



$$\mu(E) = a_{photo} \cdot \frac{1}{E^3} + a_{compton} \cdot f_{KN}(E)$$




 

Photo-effect Compton-effect

KN: Klein-Nishina Funktion





Separation of tissues



Background


Dual-Energy CT



- Simultaneous
- Spatially matched
- Retrospectively available

Background

Dual Layer CT





- Image quality improvement
- Lesion detection improvement
- Lesion characterization improvement
- Radiation dose reduction
- Contrast dose reduction
- Economic benefit






UNIKLINIK
KÖLN

Image Quality



Intra-individual comparison between abdominal virtual mono-energetic spectral and conventional images using a novel spectral detector CT



J. Doerner et al., PLoS One. 2017 Aug 24;12(8):e0183759

J. Doerner et al.: Abdom Radiol (NY). 2017 Jul 4. [Epub ahead of print]

Thorax: Image Quality


Image quality evaluation of dual-layer spectral detector CT of the chest and comparison with conventional CT imaging



CI



70 keV

40 keV
(adjusted window)



→ MonoE at low keV showed exceptional image quality compared to conventional images

Improvement of Image Quality in Unenhanced Dual-Layer CT of the Head



MonoE reconstruction

- Improves image quality
- Enhances gray-white matter contrast






UNIKLINIK
KÖLN

Lesion Detection




Assessment of arterially hyper-enhancing liver lesions using virtual monoenergetic images from spectral detector CT







HCC, S/p MWA





Diagnostic assessment of pulmonary emboli




Große Hokamp, N et al., *Journal of Computer Assisted Tomography* 2018 (accepted for publication)



UNIKLINIK
KÖLN

Breast lesions

53 yo female patient, ovary cancer



Standard IMR




MonoE 40keV




UNIKLINIK
KÖLN

Breast lesions




architectural lesion




MRI




conventional



Soft tissue metastases



Iodine Overlay



Lennartz et al., ECR 2018




Lesion characterization




Pat. w/ melanoma

cystic renal mass
DD metastasis

>not suspect






63 yo patient w/ melanoma

pulm. nodule
right upper lobe

> iodine uptake 3.1 mg/ml

> metastasis likely







UNIKLINIK
KÖLN

Pleura

Pleural fibrosis: no Iodine uptake







UNIKLINIK
KÖLN

Pleura

Pleural carcinosis: Iodine uptake






Renal Stone Characterization using Low Dose Protocols

ex vivo study with 154 stones

- Normal dose: CTDI_{vol} 10mGy
- Low dose: CTDI_{vol} 2 mGy
- Semi-automatic segmentation and pixel-by-pixel post-processing

Separation of most stones is possible irrespective of imaging protocol






Calcium Suppressed Images (CSI)

POLY

A



CaSupp

B




T1

C



STIR


D



- 34 Patients with 57 vertebral fractures
- Sensitivity 87%, Specificity 99% (PPV 95%, NPV 98%, accuracy 97%)



Reduction of Radiation and Contrast Dose







Brilliance 6, FBP, 120 kV,
100 mAs, CTDI **6,9 mGy**,
80 ml Accupaque 350, **CNR 5,3**

IQon, FBP, 120 kV,
35 mAs, CTDI **3,2 mGy**,
40 ml Accupaque 350, **CNR 4,6**

IQon, monoE 40 keV, 120 kV,
35 mAs, CTDI **3,2 mGy**,
40 ml Accupaque 350, **CNR 14,1**



Economic Benefit



HU = 39
Likert-Score: 2

Iodine = 0.3 mg/ml
Likert-Score: 5



Lesion type	n
FNH	5
Adenoma	4
Hemangioma	27
Cyst	61
HCC	11
Metastasis	29
Cirrhotic nodule	2
	139

Evaluation: benign vs. malignant / most likely diagnosis

1 = no confidence / totally unclear finding



2 = very little confidence / unclear finding

3 = moderate confidence / potential diagnosis

4 = strong confidence / likely diagnosis

5 = absolute confidence / most likely diagnosis

Diagnostic Confidence and Economic Savings



Slebocki et al.: Multiparametric Spectral Detector CT for the assessment of focal liver lesions: evaluation of diagnostic confidence



Economic Analysis

- **112 lesions** w/ score of 1-4 on conventional CT requiring follow-up scanning for confirmation of diagnosis.
- Of these, **52 lesions** changed their confidence rating to 5 on the multi-parametric dataset analysis leading to **reduction in follow-up exams by 46%**.
- **Savings of Euro 89,000 to 116,600** (MRI:US ratio of 50:50 to 75:25).



- Image quality improvement ✓
- Lesion detection improvement ✓
- Lesion characterization improvement ✓
- Radiation dose reduction ✓
- Contrast dose reduction ✓
- Economic benefit ✓

-> better patient care !



UNIKLINIK
KÖLN

david.maintz@uk-koeln.de

